

88-3 (2004)	Acetaldehyde
CAS N°: 75-07-0	EINECS N°: 200-836-8
EC-LV (8 h): - Lowest European LV (8h): 37 mg/m³ Highest European LV (8h): 180 mg/m³	EC-STLV: - Lowest European STLV: 45 mg/m³ Highest European STLV: 92 mg/m³

SUMMARY OF THE METHOD

Language: French	Reference: Aldéhydes: MétroPol Fiche 001, INRS, Paris (2005).
----------------------------	--

Summary: The method is applicable to aldehyde vapours in air. The sample is collected by pumping air through a sorbent tube containing silica gel impregnated with 2,4-DNPH. After desorption with acetonitrile the sample is analysed by HPLC with UV detection on a C18 bonded silica column.

SAMPLING

Sampler type	Pumped sorbent tube
Sampling substrate	Silica gel impregnated with 2,4-DNPH (250 - 500 mg)
Recommended flow rate	0,2 - 1 l/min
Recommended sampling time	60 - 300 min
Recommended volume	60 l

TRANSPORT AND STORAGE

Description/conditions of transport and storage incl. specific issues	Samples have to be kept refrigerated and should be analysed as soon as possible. Rapid degradation is observed.
--	---

ANALYSIS

Sample preparation	Transfer the sorbent and the glass wool into a crimp vial. Add acetonitrile (1 - 10 ml), the vial is closed and shaken for several minutes.
Analytical technique	Analysis by HPLC. An external standard (acetaldehyde DNPH-derivative) is used.

METHOD EVALUATION DATA

Range studied	-
Sampling bias	-
Analytical bias	-
Method bias	-
Sampling precision	not applicable
Analytical precision	-
Method precision	-
Limit of quantification	-
Overall uncertainty (EN 482)	< 20 % (estimated; the description file gives insufficient data for calculation).
Expanded uncertainty (prEN 482)	Insufficient data for calculation.

INFORMATION IN RELATION TO THE VALIDATION	
Is the sample dissolution procedure described widely applicable?	yes
Does the sample dissolution method include wall deposits, where applicable?	not applicable
Was a test gas atmosphere used, where applicable?	Yes (dynamic test atmosphere).
How was the recovery determined?	not mentioned
Was the sampler capacity or breakthrough volume determined?	Yes, at a concentration of 234 mg/m ³ the breakthrough volume is 26 l.
Was temperature and RH considered, where appropriate?	no
EVALUATION	
Rating category	B
Rationale for rating	Up to date methodology, interlaboratory comparisons, published work, but partially validated, brief method description, insufficient validation data in the method.
Observations	<p>In case of mixtures of aldehydes, addition of 0,03 % diethylamine to the eluent improves the chromatographic resolution.</p> <p>This method has been developed for several aldehydes, before the publication of EN 482.</p> <p>More back-up data exist for this substance but are not included in the MétroPol method description file.</p>
Similar methods	DFG